IT Manuale d'uso | EN User Manual | FR Manuel de l'utilisateur
DE Benutzerhandbuch | ES Manual de usuario | PT Manual do
utilizador

IT - Sistema di terapia a compressione d'aria - Modello: FO-3001 Manuale d'uso
EN - Air Compression Therapy System - Model: FO-3001 User Manual
FR - Système de thérapie par air comprimé - Modèle : FO-3001 Manuel de l'utilisateur
DE - Luftkompressionstherapiesystem - Modell: FO-3001 Benutzerhandbuch
ES - Sistema de presoterapia - Modelo: FO-3001 Manual de usuario
PT - Sistema de terapia de compressão de ar - Modelo: FO-3001 Manual do utilizador

IT Questo manuale è valido per la FO-3001 **EN** This manual is valid for the FO-3001
FR Ce manuel s'applique au modèle FO-3001 **DE** Dieses Handbuch ist gültig für das FO-
3001 **ES** Este manual es válido para la FO-3001 **PT** Este manual é válido para o FO-3001

INDICE – CONTENTS – SOMMAIRE INHALT – CONTENIDO - CONTEÚDO

ITALIANO	4
ENGLISH	22
FRANÇAIS	40
DEUTSCH	58
ESPAÑOL	76
PORTUGUÊS	94

ITALIANO

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ENGLISH

Be sure to read this instruction manual before operation and keep properly.

This user manual is published by Foshan Hongfeng Co., Ltd.

Foshan Hongfeng Co., Ltd does not guarantee its contents and reserves the right to improve and amend it at any time without prior notice. Amendments will however be published in a new edition of this manual.

All Rights Reserved.

FO-3001 Rev.A1 © 2023, printed in 11, 2023.

Declaration of conformity:

Foshan Hongfeng Co.,Ltd declares that FO3001 complies with following standards:

EN60601-1, EN60601-1-2, EN60601-1-11, EN62304-1
ISO10993-1, ISO10993-5, ISO10993-10, EN ISO 14971

CONTENTS

1	INTRODUCTION	24
2	INTENDED USE	24
3	CONTRANINDICATIONS	24
4	SIDE EFFECTS	24
5	SAFETY INFORMATION	25
6	SPECIFICATION	27
7	PRODUCT DESCRIPTION	28
7.1	Main unit	28
7.2	LED control panel	28
8	INSTRUCTION FOR USE	30
8.1	Sleeve use instructions	30
8.2	Operation instructions	30
9	CLEANING AND MAINTENANCE	32
9.1	Cleaning and care	32
9.2	Maintenance	32
10	TROUBLESHOOTS	33
11	STORAGE	33
12	DISPOSAL	33
13	SYMBOLS	33
14	INFORMATION ON ELECTROMAGNETIC COMPATIBILITY (EMC)	35

1 INTRODUCTION

Air Compression Therapy System consists of air pressure sensor, air pump, sleeves etc working together as one unit. The air pump is connected to the dedicated sleeves via a series of hoses. The compression massage direction is from limb end to body center by inflating the air chambers sequentially and then deflating as one cycle, the pressure can be adjusted to avoid any discomfort to the patient. The sleeve works under the action of an analogue and mechanical regulation.

Air Compression Therapy System, in medical market, it has a sequential squeezing from distal to proximal, thus help to improve the circulation of blood and lymph, and to prevent the DVT and relieve lymph edema.

2 INTENDED USE

The product is indicated for use by medical professionals and patient at home, who are under medical supervision, reducing swelling and preventing thrombosis in the lower extremities or treatment of truncal or arm breast cancer-related lymphedema. such as: Primary lymphedema, edema following trauma and sport injures, post-immobilization edema, Venous insufficiencies, Lymphedema.

The system can be used at home or in hospital, to patients who are over 18 years old.

3 CONTRANINDICATIONS

- (1) Acute inflammatory skin diseases
- (2) Arrhythmia
- (3) Erysipelas
- (4) Deep thrombophlebitis
- (5) Pulmonary edema
- (6) Acute venous thrombosis
- (7) Unstable hypertension
- (8) The installation of artificial pacemaker patients.
- (9) Presence of the arterial insufficiency

4 SIDE EFFECTS


- Feel discomfort when the output pressure is too large.
- Redness, itching, or discomfort
- Discomfort and pain sensations during treatment with a pressure of 120mmHg



5 SAFETY INFORMATION

It is important that you read all the warning and precautions included in this manual because they are intended to keep you safe, prevent injury and avoid a situation that could result in damage to the device.

- Do not use this device if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death.
- Not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- The devices are not intended for unattended personal use by patients who have noncompliant, emotionally disturbed, dementia, or low IQ.

Repair is not allowed. When you found the machine abnormal, please do not disassemble it casually or privately. 	
If feel abnormal or painful when using, please stop using the product immediately.	For the first time, the treatment time should not exceed 20 minutes. For the first time, you should consult your doctor and treatment less than 20 minutes. Then gradually extend the treatment time according to your doctor's suggestions (long-term use may have side effects).
Store the equipment away from infants and young children. May cause harm such as electric shock.	Be careful! Keep water or foreign matter out of the machine.
The equipment should not be placed in places where it is difficult to disconnect the power supply.	Only use 220-240V supply connection, please unplug the power correctly after finished treatment. Inappropriate supply connection may cause device damage.

Keep the hand dry when plugging in or unplugging the power supply. Wet hand may cause electric shock 	Power supply cord cannot bear the excessive stress, cannot be lain under the tables or chairs. Broken may cause electric shock. 
Do not bend and fold the power supply cord. 	Do not make the garments near the heat of the metal, stoves, cigarettes and so on. This will cause damage to garments. 
Do not use oil, benzene, alcohol, gasoline and chemical agents to wash the garments. Otherwise, the lifetime of the garments will be shortened. Please wipe garments with damp cloth.	Do not unzip the garment when the unit is working on. This may cause damage to the garment. Always turn off the power before unzipping the garments.
Do not use the machine for nudity. Wear underwear when using. 	

6 SPECIFICATION

ITEM	TECHNICAL INDEX
Number of chambers	4 chambers
Machine size	240 × 200 × 120mm
Weight (with accessories)	2.22kg
Control way	Manual
Treatment time	10~30 minutes
Treatment pressure	30~250mmHg
Pressure error	<15%
Cycle time	60s, ± 1s
Treatment modes	1 mode
Chamber Setting	No Settable
Input voltage	AC 220-240V, 50-60Hz
Power consumption	65VA
Noise level	≤55dB
Protection against electric shock	Class II, Type BF Applied Part
Mode of operation	Continuous operation
Degree of waterproof	IP21
Product life	5 years
Operating environment	+10°C - +40°C, 10% - 95%RH 86kPa-106kPa
Storage environment	-40°C - +70°C, 10% -100%RH 50kPa-106kPa
Sleeves size	Leg cuff: 110*70cm Arm cuff: 90*60cm Waist cuff: 130*37cm
Extension part	Width 10cm (for leg sleeve)


7 PRODUCT DESCRIPTION

7.1 Main unit

The product is mainly composed of machine, sleeves and connection hoses etc.

7.2 Control panel


- Pressure: pressure can be adjusted between 30 mmHg and 250 mm Hg (pressure error +/- 15%).
 - To increase pressure, rotate the knob clockwise;
 - To reduce pressure, rotate the knob anticlockwise;
 - Perform 2 inflating cycles to minimum and then adjust pressure.
- Power: between the pressure knob and the time knob, there is a green (0 - I) OFF-ON button to start the device.
- Time: time can be adjusted between 0 and 30 minutes. After 30 minutes of machine operation, it is necessary to wait at least 30 minutes (with the machine off) before starting a new session.

- Caution: if you turn the knob to “Always On” automatic power off is disabled. Therefore, the machine must be switched off manually.
- To increase time, rotate the knob clockwise.

7.2.1 Mode description

Sequential mode: pressure is exerted by only one chamber at a time; this means that the 1st chamber is inflated and after the latter has deflated, the inflation and then deflation of the 2nd begins and so on up to the 4th air chamber and then starts again from the 1st.

In the case of the leg, inflation starts from the air chamber at the foot, in the case of the cuff it starts from the air chamber at the hand and in the case of the band it starts from the air chamber at the upper part of the abdomen .

8 INSTRUCTION FOR USE

8.1 Sleeve use instructions

- i. Insert the power plug.
- ii. Connect the hose to sleeve(s) and put it (them) on body.

- Leg Garment

In order to avoid knee pain, please don't fold sleeve when use.

- Wear thin clothes and relax naturally.
 - Pull the zipper up to the top.
 - Face to control panel of the machine, make the air outlets of sleeves upwards, and connect the connector of hose to the air outlet of machine.
- iii. Turn on the machine, set pressure level and then the treatment time.

8.2 Operation instructions

BEFORE CONNECTING THE POWER CORD AND CARRYING OUT ANY OPERATION ON THE MACHINE, IT IS ABSOLUTELY NECESSARY TO REMOVE THE SCREW THAT BLOCKS THE MOTOR OF THE DEVICE.

Note: After receiving the machine, please open filter cover, and remove the security screw before use it. The screw is **ONLY** used for fixing and protecting the machine from damage during transportation.

After this important operation:

- 1) Please find a flat plane to place the machine.
- 2) For installation, it just needs to connect the host connector to the machine and sleeve, then plug in supply mains.

Step 1 Insert the main plug of power supply cord into electrical outlet of 220~240V.

Step 2 Check if the hose has been connected well to the sleeve(s): if yes, skip to step 3; if not, please fit the nozzle of the hose into the outlets sleeve(s).

Step 3 Connect the main end of the connection hose to the outlet of main unit.

Step 4 Check and ensure all connections to the main unit and sleeve firmly.

Step 5 Put on the sleeves.

CAUTION:


The connection hose cannot place on the neck to avoid strangulation or injure the neck.

Step 6 Set the power button on “I” position to turn on the machine.


Note: Set the power button on “O” position to turn off.

Step 7 Set the pressure by rotating the knob clockwise to increase and anticlockwise to reduce.

Note: for first time use, it had better start with the pressure not too high, maybe 120 or 150 mmHg, and then could adjust the pressure to the level that makes you feel comfortable.

Step 8 Set the time by turning the knob clockwise to start inflation and increase the time.

Note: time range is 10~30 minutes. if you turn the knob to “Always On” automatic power off is disabled. Therefore, the machine must be switched off manually.

CAUTION:

Do not unzip the sleeve(s) during the treatment. This may cause damage the sleeve(s). Always turn off the device before unzipping.

Had better not disconnect hose from the sleeve frequently to ensure the connection tight.

9 CLEANING AND MAINTENANCE

Fully comply with the following daily maintenance requirements are necessary, to make sure your product operate intact and guarantee the product long-term performance and safety.

9.1 Cleaning and care

Not to clean the device during treatment, be sure that the device is turned off before cleaning the device.

Clean the device with a soft, slightly damp cloth. In case of heavier dirt build-up, you may also add a mild detergent.

Ensure that no water enter into the machine. If this happens, only use the device after it is completely dry.

Do not subject the device to moisture or dampness. Do not hold the device under running water, do not submerge it in water or other liquids.

For reasons of hygiene, each user should use his/her own set of sleeves. Do not use any organic solvent or aggressive reagent for sleeve cleaning. Do not wash sleeve with water.

Carefully clean the surface of the sleeves with a damp cloth.

9.2 Maintenance

- 1) The manufacturer didn't authorize any maintenance agencies abroad. If your device has any problem, please contact the distributor. The manufacturer will not be responsible for the results of maintenance or repairs by unauthorized persons.
- 2) The user must not attempt any repairs to the device or any of its accessories. Please contact the distributor for repair.
- 3) Opening of the equipment by unauthorized agencies is not allowed and will terminate any claim to warranty.

Each product in manufacturing has been passed systematic validation. The performance is stable and does not need to undertake calibration and validation.

If your product can't reach the expected performance and the basic function has changed in normal use, please contact the distributor.

10 TROUBLESHOOT

Common Malfunctions	Check and Confirm
If the machine cannot start	♦ if the plug is connected well to power ♦ if the machine is turned on ♦ if the power supply is normal (220-240V) ♦ if the machine automatically turned off when treatment time is over.
There are leaks and other strange sounds	♦ if the connection pipe and connector are damaged ♦ if the connector is connected to the machine well ♦ if the connecting tube is pressed or bent

11 STORAGE

Do not make any sharp kinks in the connection hoses or sleeves. After use, put the sleeves in original package. Do not expose the device to direct sunlight, protect it against dirt and moisture. Store the machine in a cool, well-ventilated place. Never place any heavy objects on the machine.

12 DISPOSAL

At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health. You are obligated to dispose of device correctly. Consult your municipal authority or your dealer for information about disposal.

13 SYMBOLS

Electrical devices are recyclable material and should not be disposed of with household waste after their useful life! Help us to protect the environment and save resources and take this device to the appropriate collection points. Please contact the organization which is responsible for waste disposal in your area if you have any questions.

Type of protection against electric shock: Class II Equipment

Degree of protection against electric shock: Type BF Applied Part

Refer to User Manual

IP21

The first number **2**: Protected against solid foreign objects of $\Phi 12.5\text{mm}$ and greater. The second number **1**: Protected against vertically falling water drops when vertically dripping. Vertically falling drops shall have no harmful effects when the enclosure is tilted at vertically dripping, on either side of the vertical.

Serial number

Year			Month		PO		Serial number		

Manufacturer information


EC Authorized Representative

Manufacture date

CE Marking, Certificate issued by NB No.:0598

Up to 4 layers

Fragile items, be careful

Product fear of rain, keep dry

14 INFORMATION ON ELECTROMAGNETIC COMPATIBILITY (EMC)

- 1) This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information below, and this product can be affected by portable and mobile RF communications equipment.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields near the product. This may result in incorrect operation of the product.
- 3) This product has been thoroughly tested and inspected to assure proper performance and operation!
- 4) This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used

Guidance and manufacture's declaration – electromagnetic emission		
FO3001 is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	FO3001 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacture's declaration – electromagnetic immunity			
FO3001 is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor is covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	±1kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the DEVICE requires continued operation during power mains interruptions, it is recommended that the DEVICE be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	10V/m	10V/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacture's declaration – electromagnetic immunity			
FO3001 is intended for use in the electromagnetic environment specified below. The user should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	10V/m & table 9	10V/m & table 9	Portable and mobile RF communications equipment should be used no closer to any part of the Device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.167\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.333\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 3) At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 4) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. ^c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Device is used exceeds the applicable RF compliance level above, the Device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Device. ^d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [Vi] V/m.			

Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment (Table 9)						
Test freq. (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	Immunity Test Level (V/m)
385	380-390	TETRA 400	Pulse modulation ^{b)} 18Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM ^{c)} ±5kHz deviation 1kHz sine	2	0.3	28
710	704-787	LTE Band 13, 17	Pulse modulation ^{b)} 217Hz	0.2	0.3	9
745						
780						
810	800-960	GSM800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18Hz	2	0.3	28
870						
930						
1720	1700-1990	GSM1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation ^{b)} 217Hz	2	0.3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217Hz	0.2	0.3	9
5500						
5785						
If necessary to achieve the immunity test level, the distance between the transmitting antenna and the FO3001 may be reduced to 1m. The 1m test distance is permitted by IEC 61000-4-3.						
d	For some services, only the uplink frequencies are included.					
e	The carrier shall be modulated using a 50 % duty cycle square wave signal.					
f	As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.					

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MedNet EC-REP GmbH